

# Appendix 7: Decontamination of reusable non-invasive care equipment



**Infection prevention and control team/health protection team:**

Name: \_\_\_\_\_

Designation: \_\_\_\_\_


Contact Number: \_\_\_\_\_

Routine decontamination of reusable non-invasive care equipment.

- Check manufacturer's instructions for suitability of cleaning products especially when dealing with electronic equipment.
- Wear appropriate PPE eg disposable, non-sterile gloves and apron.

Is equipment contaminated with urine/vomit/faeces or been used on a patient with a known/suspected infection/colonisation?

Is equipment contaminated with blood or body fluid as in box 1?



- Follow manufacturer's instructions for safe preparation, dilution, and contact time.
- Clean the piece of equipment from the top or furthest away point.
- Discard disposable cloths/paper roll immediately into the healthcare waste receptacle.
- Discard detergent/disinfectant solution in the designated area.
- Clean, dry and store re-usable decontamination equipment.
- Remove and discard PPE.
- Perform hand hygiene.

- BOX 1**
- Cerebrospinal fluid
  - Peritoneal fluid
  - Pleural fluid
  - Synovial fluid
  - Amniotic fluid
  - Semen
  - Vaginal secretions
  - Breast milk
  - Any other body fluid with visible blood (excluding urine)