

Trauma  
 Rehabilitation and Disability  
 Complex Disability Equipment - Prosthetic Specialised Services For People Of All Ages With Limb Loss

D01/S/D

Code	Metric	Numerator	Denominator	Period Type	Frequency	Data Source Numerator	Data Source Denominator	Target	Notes	Q1	Q2	Q3	Q4
PROS01	Proportion of clinically appropriate patients who are cast for primary prosthesis within 2 weeks of their primary appointment	Of those in the denominator, the number of patients who are cast for primary prosthesis within 2 weeks of the appointment	The total number of clinically appropriate patients having their primary appointment in the reporting period	Quarter	Quarterly	Provider	Provider	80	<p>Clinically appropriate is defined as the patient is deemed by the MDT (2 or more clinicians) to be a suitable candidate for prosthetic provision, and is physically and psychologically ready for casting and proceeding to rehabilitation.</p> <p>If a patient is deemed not clinically appropriate for any prosthetic provision at their primary / first appointment, and this is unlikely to change, the reasons for this decision must be explained to the patient and relatives / carers, and access to the service as a non limb user offered.</p> <p>If suitability is not clear, or the patient is suitable but not yet ready at their primary / first appointment with the prosthetic centre post amputation, and then the clock can be stopped for the 2 week timescale for casting, and restarted if indicated at the point when the MDT agrees it is clinically appropriate to proceed.</p> <p>If the clock is stopped, a suitable timepoint for review of suitability for prosthesis provision must be put in place.</p>	Apr-22 - Jun-22	Jul-22 - Sep-22	Oct-22 - Dec-22	Jan-23 - Mar-23
PROS02	Proportion of patients who have a recorded Socket Comfort Score of >7 at each fitting	Of those in the denominator, the number of patients who have a Score of >7 recorded	The total number of patients having a Socket Comfort Score recorded in the reporting period	Quarter	Quarterly	Provider	Provider	80	<p>Socket comfort score is a validated numerical rating scale, where 0 represents the lowest level of comfort (least comfortable) and 10 the highest level of comfort (most comfortable) related to comfort of prosthetic socket fit.</p> <p>Count patients more than once if they have had more than 1 fitting in the reporting period.</p>	Apr-22 - Jun-22	Jul-22 - Sep-22	Oct-22 - Dec-22	Jan-23 - Mar-23
PROS03	Proportion of new sockets delivered within 4 weeks of the initial cast	Of those in the denominator, the number which are delivered within 4 weeks of the initial cast taking place	The total number of new sockets being initially cast in the reporting period	Quarter	Quarterly	Provider	Provider	80	<p>If patient is not able for clinical or other reasons to proceed with appointments for progression of the socket, then the clock can be stopped and restarted when appropriate.</p>	Apr-22 - Jun-22	Jul-22 - Sep-22	Oct-22 - Dec-22	Jan-23 - Mar-23
PROS04	Proportion of new sockets delivered within 4 appointments (to include the initial cast appointment and the delivery appointment)	Of those in the denominator, the number delivered within 4 appointments	The total number of new sockets being delivered in the reporting period	Quarter	Quarterly	Provider	Provider	80	<p>If patient is not able for clinical or other reasons to proceed with appointments for progression of the socket, then the clock can be stopped and restarted when appropriate.</p>	Apr-22 - Jun-22	Jul-22 - Sep-22	Oct-22 - Dec-22	Jan-23 - Mar-23
PROS05	Proportion of active patients who are offered an annual review (may be carried out by telephone if appropriate)	Of those in the denominator, the number who are offered a review	The total number of active patients who are due an annual review in the reporting period	Quarter	Quarterly	Provider	Provider	90	<p>An active patient is one that has had input from the service at least once in the last 2 years.</p>	Apr-22 - Jun-22	Jul-22 - Sep-22	Oct-22 - Dec-22	Jan-23 - Mar-23